

CANES & CRUTCHES

CANES



CANES WITH HEIGHT ADJUSTMENT FASHION

Aluminum canes facilitate walking for elderly or disabled persons. They are equipped with height adjustment function and with handgrip made of wood. Canes are available in 2 colors: black & brown.



Height adjustment (range)	Height adjustment steps	Weight	Internal diameter of rubber foot	Diameter of rubber foot's base	Max. weight of user
71-94 cm	2,5 cm	0,4 kg	19 mm	38 mm	do 100 kg



CANES WITH HEIGHT ADJUSTMENT

Aluminum, fashion canes facilitate walking for elderly or disabled persons. They are equipped with height adjustment function and with handgrip made of maple wood. Canes are available in 4 styles: scotch, flowers, chessboard & graphite (marbled).



Height adjustment (range)	Height adjustment steps	Weight	Internal diameter of rubber foot	Diameter of rubber foot's base	Max. weight of user
71-94 cm	2,5 cm	0,4 kg	19 mm	38 mm	do 100 kg