

# WALKERS & ROLLATORS

## WALKERS



RECIPROCATING  
FOLDABLE WALKER



VERSION WITH TWO WHEELS

The reciprocating function may be blocked and the walker can be used as a rigid one.



## RECIPROCATING REHABILITATION WALKERS

Rehabilitation walkers are mobility aids designed to support walking and standing for adult and junior patients with mobility problems. Walkers can be used for walking self-reeducation after disease or injury. Height of a walker can be adjusted to fit each user's height. Ergonomic and anatomically shaped handgrips deliver proper stabilization.

Frame material	Width	Length	Push handles height (range)	Height adjustment steps	Max. weight of user	Weight
aluminum	62 cm	53,5 cm	82,3-100 cm	2,5 cm	136 kg	3 kg