

STAIRCLIMBERS & RAMPS FOR THE DISABLED

RAMPS



TR/LR PORTABLE WHEELCHAIR RAMPS (telescopic)

This type of telescopic ramps is very easy to open and to fold. These are very universal products used for negotiating curbs, but also for rolling heavy loads in to cars. It is possible to adjust their spacing freely for wheelchairs or scooters with different widths. Wide variety of available sizes allows to climb up different obstacles, from the lowest curbs up to 6 steps.



	TR-150	TR-210	TR-300
Width x Length	20 x 150 cm	20 x 210 cm	20 x 300 cm
Max. load	225 kg	225 kg	225 kg
Weight	2 x 5,5 kg	2 x 7,5 kg	2 x 9,5 kg



SC PORTABLE WHEELCHAIR RAMPS (with ribbed surface)

Ribbed surface ramps are new version of our well known foldable ramps. They have new rolling surface, made of ribbed aluminium. Thanks to that they offer even better traction and (specially) easy maintenance and cleaning. They are single folding cases available in 2 sizes.



	SC-120	SC-150
Width x Length	71 x 120 cm	71 x 150 cm
Max. load	300 kg	300 kg
Weight	9,1 kg	10,9 kg